


## Updating Microphones Guide (PureTest)

While using your PureTest app, you may notice a label next to your microphone indicating that it is about to expire, or has already expired. If you notice this, you will need to contact SHOEBOX support to have your microphone recalibrated and sent back to you. It's important to be proactive in contacting SHOEBOX for microphone recalibrations, as using expired hardware may affect both your compliance with applicable regulations and the accuracy of your test results.

Once you have received your newly calibrated microphone from SHOEBOX, you will need to add it in the app to replace your old, expired microphone.

### To update your microphone:

1. From the home screen, tap "Settings" in the bottom right corner to access the Settings screen
2. Under the Hardware heading, tap on, "Headphones and Microphone"  
The Microphone panel will display the name of the microphone that is currently installed on your iPad, and will indicate if the microphone is expired or about to expire
3. Tap **Remove** in the top right of the microphone panel to remove the old microphone
4. Tap on the **+** icon to open the "Change Microphone" panel
5. Install your newly calibrated microphone using one of the following options:
  - a) Select your new microphone from the bulleted list of available microphones
  - b) Enter the 6-digit ID number etched on the side of your new microphone



**Once installed, your new microphone will be listed on the "Microphone" panel**