

SHOEBOX PureTest vs. SHOEBOX Standard/Pro: Key Differences

| Feature | PureTest | SB Standard/Pro |
|--|---|---|
| Fleet Management - settings locked from Portal | ✓ All test settings and workflow steps controlled from the web portal | ✗ All test settings are set in the app |
| Enforce checks | ✓ Functional check and Room Scan can be enforced to block testing | ✗ Only the room scan can be enforced |
| Multi-user Login | ✓ Each user may have a login under a single license | ✗ Not available, only one user may be logged in at a time |
| Offline Shift Analysis | ✓ Fully offline shift analysis available by selecting Projects for offline testing. | ✗ Not available, must be online. |
| Workflow | ✓ Wizard takes you through the entire testing process | ✗ Steps must be selected to proceed through workflow |
| Single Results Document | ✓ Results, baseline, shift analysis, hearing grade overlay and signatures can be viewed immediately after testing and does not require a data connection. | ✗ A basic audiogram can be viewed immediately after testing. To view the summary report, it must be downloaded using a data connection. |
| Data Privacy | ✓ Results are stored hidden until synced and then removed from the cache. | ✓ Audiograms stay on the device until manually deleted. They are protected behind user logins and passcodes. |
| Languages | ✓ Fully localized in English, French, and Spanish | ✓ Available in 27 subject-facing languages, and 4 Examiner-facing languages |